

# Primer Designing for AKT1 gene using Primer-BLAST

## 1. Get the Exon and Intronic Sequences from NCBI

- **Search AKT1 Gene:** Find the **AKT1 gene** sequence on **NCBI** (e.g., **RefSeq ID: NG\_012188.1**).
  - **Select Exon Sequence:** Choose the specific exon region you want to design primers for.
  - **Add Introns:** Include introns at the **5' and 3' ends** to ensure accurate amplification and avoid mismatch.
- 

## 2. Input Sequence into Primer-BLAST

- **Go to Primer-BLAST:** Access the [Primer-BLAST Tool](#).
  - **Enter Sequence:** Paste the full sequence (including exon and intron regions) into the "Enter target sequence" box.
  - **Select Organism:** Choose **Homo sapiens** for the AKT1 gene.
- 

## 3. Set Primer Parameters in Primer-BLAST

- **Primer Size:** Set primer length to **18-25 bp**.
  - **Amplicon Size:** Set the amplicon size to **100-300 bp** for standard PCR.
  - **Melting Temperature (T<sub>m</sub>):** Set T<sub>m</sub> to **58-60°C** to optimize primer binding.
  - **GC Content:** Choose primers with **40-60% GC content**.
- 

## 4. Primer Specificity Settings

- **Primer Database:** Select **nr** (Non-Redundant) or **nt** (Nucleotide) database.
  - **Product Size:** Ensure the product size falls within the desired range (**100-300 bp**).
  - **Check for Non-Specific Binding:** Ensure the "Allow Primer3 to evaluate the potential for non-specific amplification" option is checked.
-

## 5. Generate Primers

- Click **"Get Primers"**: Primer-BLAST will generate forward and reverse primer pairs based on your sequence.
- 

## 6. Result Interpretation

Once Primer-BLAST generates the primer pairs, here's how to interpret the results:

- **Primer Sequences**: You will receive both **forward** and **reverse** primer sequences. Ensure they match your desired gene region, and that the **amplicon size** falls within the expected range (e.g., **100-300 bp**).
- **Amplicon Size**: Check the **amplicon size** (product size) reported in the result. Ensure it matches your expected target (e.g., **150-250 bp** for efficient PCR amplification).
- **T<sub>m</sub> (Melting Temperature)**: Primer-BLAST will display the **T<sub>m</sub>** values for both primers. The T<sub>m</sub> of both primers should be within **±2°C** of each other for efficient binding during PCR. For example, if your forward primer has a T<sub>m</sub> of **59°C**, the reverse primer should be between **57°C and 61°C**.
- **GC Content**: Ensure the **GC content** of both primers is within the optimal range of **40-60%**. This helps in creating stable primer-template bindings.
- **Specificity Report**:
  - **No Non-Specific Binding**: Ensure the primers are specific to your target region and do not bind to other genomic sequences. Primer-BLAST will provide a **specificity check** and show if any primer binds to unintended regions. If non-specific binding is detected, you may need to adjust the primer design (e.g., extend or reduce the region of amplification).
- **Secondary Structures**: Primer-BLAST will flag any **hairpins** or **primer dimers**. A **hairpin** is a secondary structure where a primer binds to itself, while **primer dimers** occur when the two primers bind to each other. Both can reduce PCR efficiency. Ensure neither structure is present for both forward and reverse primers.
- **Self-Dimer and Cross-Dimer Check**: Ensure that primers do not form **self-dimers** or **cross-dimers** (binding with each other). If dimerization occurs, adjust the primer sequences by changing the primer length, GC content, or sequence composition.